

Anti-Melan-A/Mart-1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q16655

Item Number: CM0024868



Introduction

Rabbit polyclonal antibody targeting Melan-A/MART-1, a melanosome biogenesis protein. Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for melanoma and melanocyte studies.

[Learn More](#)

Parameter	Value
Product Name	Melan-A/MART-1 Rabbit pAb
Remarks/Alias	MART1; MART-1; Melan-A/MART-1
Species	Human
GeneID (human)	2315
GeneID	2315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 48-118 of human Melan-A/MART-1 (NP_005502.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16655
Protein Weight	13kDa
Shipping	Ice bag