

Anti-Melan-A/Mart-1 Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - Q16655

Item Number: CM0007214



Introduction

Melan-A/MART-1 Rabbit mAb (Q16655) validated for WB, IF, IHC, ELISA. Specific for human melanoma antigen, 13 kDa. Ideal for melanocyte and melanoma research.

[Learn More](#)

Parameter	Value
Product Name	Melan-A/MART-1 Rabbit mAb
Remarks/Alias	MART1; MART-1; Melan-A/MART-1
Species	Human
GenelD (human)	2315
GenelD	2315
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 11-110 of human MLANA (NP_005502.1)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q16655
Protein Weight	13kDa
Shipping	Ice bag