


Anti-Melk Rabbit Polyclonal Antibody For Wb And Elisa - Q14680

Item Number: CM0013597



Introduction

Anti-MELK rabbit polyclonal antibody for Western blot and ELISA applications in human, mouse, and rat samples. Targets maternal embryonic leucine zipper kinase involved in cell cycle and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	MELK Rabbit pAb
Remarks/Alias	HPK38; MELK
Species	Human
GeneID (human)	9833
GeneID	9833
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 412-651 of human MELK (NP_055606.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14680
Protein Weight	75kDa
Shipping	Ice bag