


Anti-Melk/Hpk38 Rabbit Monoclonal Antibody For Wb, Elisa - Q14680

Item Number: CM0005400



Introduction

High-quality rabbit monoclonal antibody against human MELK (HPK38). Validated for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for cell cycle, apoptosis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MELK Rabbit mAb
Remarks/Alias	HPK38; MELK
Species	Human
GeneID (Human)	9833
GeneID	9833
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MELK (NP_055606.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14680
Protein Weight	75kDa
Shipping	Ice bag