


Anti-Memo1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y316

Item Number: CM0030717



Introduction

Rabbit polyclonal antibody targeting human MEMO1 for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studying ERBB2-driven cell motility and breast cancer migration.

[Learn More](#)

Parameter	Value
Product Name	MEMO1 Rabbit pAb
Remarks/Alias	MEMO; C2orf4; CGI-27; NS5ATP7; MEMO1
Species	Human
GeneID (human)	51072
GeneID	51072
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-297 of human MEMO1 (NP_057039.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y316
Protein Weight	34kDa
Shipping	Ice bag