


Anti-Menin Polyclonal Antibody For Wb, Elisa - 000255

Item Number: CM0013063



Introduction

Rabbit polyclonal antibody targeting human Menin, validated for WB and ELISA. Cross-reacts with mouse. Menin is a key transcriptional regulator and component of the MLL/SET1 histone methyltransferase complex.

[Learn More](#)

Parameter	Value
Product Name	Menin Rabbit pAb
Remarks/Alias	MEAI; SCG2; Menin
Species	Human
GeneID (human)	4221
GeneID	4221
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 266-615 of human Menin (NP_000235.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00255
Protein Weight	68kDa
Shipping	Ice bag