

Anti-Mep1A Polyclonal Antibody For Wb And Elisa - Q16819

Item Number: CM0030326



Introduction

MEP1A Rabbit pAb is a rabbit polyclonal antibody against human meprin A subunit alpha. Validated for WB and ELISA, it detects MEP1A in human, mouse, and rat samples. Ideal for metalloprotease research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	MEP1A Rabbit pAb
Remarks/Alias	PPHA; MEP1A
Species	Human
Gene ID (Human)	4224
Gene ID	4224
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 440-700 of human MEP1A (NP_005579.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16819
Protein Weight	84kDa
Shipping	Ice bag