

Anti-Mepe Polyclonal Antibody For Wb, Elisa - Q9Nq76

Item Number: CM0026482

Introduction

High-quality rabbit polyclonal antibody against human MEPE, suitable for WB and ELISA. Detects MEPE in human, mouse, and rat samples. Ideal for bone mineralization and dental pulp research.

[Learn More](#)

Parameter	Value
Product Name	MEPE Rabbit pAb
Remarks/Alias	OF45; MEPE
Species	Human
GeneID (Human)	56955
GeneID	56955
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 90-256 of human MEPE (NP_064588.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQ76
Protein Weight	58kDa
Shipping	Ice bag