

Anti-Mers-Cov Nucleoprotein Polyclonal Antibody For Wb, Elisa - K9N4V7

Item Number: CM0019859



Introduction

Rabbit polyclonal antibody against MERS-CoV Nucleoprotein (N protein), validated for WB and ELISA. Targets the nucleocapsid protein essential for viral RNA packaging and replication. Suitable for coronavirus research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	MERS-CoV Nucleoprotein Rabbit pAb
GeneID	14254601
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-413 of coronavirus Nucleoprotein (YP_009047211.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	MERS-CoV
SWISS	K9N4V7
Protein Weight	45kDa
Shipping	Ice bag