

Anti-Mers-Cov Nucleoprotein Rabbit Monoclonal Antibody For Wb, Elisa - K9N4V7

Item Number: CM0006593



Introduction

Rabbit monoclonal antibody targeting MERS-CoV nucleoprotein, validated for WB and ELISA applications. Recognizes MERS-CoV specifically, with a predicted molecular weight of 45 kDa. Ideal for coronavirus research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	MERS-CoV Nucleoprotein Rabbit mAb
Species	Human
Gene ID (Human)	14254601
Gene ID	14254601
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-413 of coronavirus Nucleoprotein (YP_009047211.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	MERS-CoV
SWISS	K9N4V7
Protein Weight	45kDa
Shipping	Ice bag