


Anti-Mers-Cov Spike Rbd Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - K9N5Q8

Item Number: CM0024482



Introduction

Rabbit polyclonal antibody targeting MERS-CoV Spike RBD (YP_009047204.1). Validated for WB, IF/ICC, ELISA; reacts with MERS-CoV. Suitable for viral entry and spike protein research.

[Learn More](#)

Parameter	Value
Product Name	MERS-CoV Spike RBD Rabbit pAb
Species	Coronavirus
GeneID	14254594
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-725 of coronavirus Spike RBD (YP_009047204.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	MERS-CoV
SWISS	K9N5Q8
Protein Weight	149kDa
Shipping	Ice bag