


Anti-Mertk Monoclonal Antibody For Wb - Q12866

Item Number: CM0002925



Introduction

High-quality rabbit monoclonal antibody against human MERTK, validated in WB and ELISA. Targets the receptor tyrosine kinase Mer, essential in macrophage efferocytosis and retinal health. SWISS: Q12866.

[Learn More](#)

Parameter	Value
Product Name	MERTK Rabbit mAb
Remarks/Alias	MER; RP38; c-Eyk; c-mer; Tyro12; MERTK
Species	Human
Gene ID (Human)	10461
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MERTK (Q12866)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q12866
Protein Weight	110kDa
Shipping	Ice bag