


Anti-Mertk/Mer Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q12866

Item Number: CM0028130



Introduction

High-quality rabbit polyclonal antibody against human MERTK (MER, Tyro12), validated for WB, IHC-P, IF/ICC, ELISA. Detects MERTK in human, mouse, and rat tissues; ideal for phagocytosis and cancer signaling studies.

[Learn More](#)

Parameter	Value
Product Name	MERTK Rabbit pAb
Remarks/Alias	MER; RP38; c-Eyk; c-mer; Tyro12; MERTK
Species	Human
Gene ID	10461
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-499 of human MERTK (NP_006334.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12866
Protein Weight	110kDa
Shipping	Ice bag