

## Anti-Mesdc2 Polyclonal Antibody For Wb And Elisa - Q14696

Item Number: CM0018278

### Introduction

Rabbit polyclonal antibody against MESDC2 (LRP chaperone MESD, OI20), validated for WB and ELISA. Recognizes human and mouse MESDC2. Recombinant immunogen (aa1-234). Ideal for LDLR and Wnt signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | MESDC2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | BOCA; OI20; MESDC2                                                                                                |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 23184                                                                                                             |
| GeneID           | 23184                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human MESDC2 (NP_055969.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | Q14696                                                                                                            |
| Protein Weight   | 26kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |