


Anti-Mesothelin Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13421

Item Number: CM0005286



Introduction

Recombinant rabbit monoclonal antibody against human Mesothelin (MSLN/MPF), validated for Western blot, IHC-P, and ELISA. Targets amino acids 296-606. Suitable for cancer and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Mesothelin Rabbit mAb
Remarks/Alias	MPF; SMRP; Mesothelin
Species	Human
Gene ID (Human)	10232
Gene ID	10232
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 296-606 of human Mesothelin (NP_037536.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q13421
Protein Weight	67KDa/68KDa/71KDa
Shipping	Ice bag