


Anti-Mesothelin Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q13421

Item Number: CM0013496



Introduction

Rabbit polyclonal antibody against human Mesothelin (MSLN, MPF, SMRP) validated for WB, IF-P, IHC-P, ELISA. Reacts with Human, Mouse, Rat. Suitable for cancer and mesothelial research.

[Learn More](#)

Parameter	Value
Product Name	Mesothelin Rabbit pAb
Remarks/Alias	MPF, SMRP, Mesothelin
Species	Human
GeneID (Human)	10232
GeneID	10232
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 296-606 of human Mesothelin (NP_037536.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13421
Protein Weight	69kDa
Shipping	Ice bag