


Anti-Mesothelin Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q13421

Item Number: CM0007824



Introduction

High-quality Anti-Mesothelin Rabbit Monoclonal Antibody validated for WB, IF-P, ELISA. Reacts with human and rat. Ideal for mesothelin/cancer research. This antibody detects mesothelin, a cell adhesion protein, and is suitable for diagnostic assay development. SWISS: Q13421.

[Learn More](#)

Parameter	Value
Product Name	Mesothelin Rabbit mAb
Remarks/Alias	MPF; SMRP; Mesothelin
Species	Human
Gene ID	10232
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 500-600 of human Mesothelin (Q13421)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13421
Protein Weight	69kDa
Shipping	Ice bag