


Anti-Mesothelin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13421

Item Number: CM0002561



Introduction

Rabbit monoclonal antibody against human Mesothelin (MSLN), validated for WB, IHC-P, and ELISA. Specific for the CAK1/MPF tumor biomarker, ideal for cancer research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Mesothelin Rabbit mAb
Remarks/Alias	MPF; SMRP; Mesothelin
Species	Human
GeneID (Human)	10232
GeneID	10232
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 296-606 of human Mesothelin (NP_037536.2).
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q13421
Protein Weight	69kDa
Shipping	Ice bag