


Anti-Mesothelin Rabbit Polyclonal Antibody For Wb, Elisa - Q13421

Item Number: CM0027954



Introduction

Rabbit polyclonal antibody against human Mesothelin (SWISS: Q13421) for Western blot and ELISA applications. Recognizes human and mouse Mesothelin; useful for cancer and megakaryocyte research.

[Learn More](#)

Parameter	Value
Product Name	Mesothelin Rabbit pAb
Remarks/Alias	MPF; SMRP; Mesothelin
Species	Human
GeneID	10232
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 296-606 of human Mesothelin (NP_037536.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13421
Protein Weight	69kDa
Shipping	Ice bag