

Anti-Mest Rabbit Polyclonal Antibody For Wb, If/Icc, And Elisa - Q5Eb52

Item Number: CM0025254

Introduction

Rabbit polyclonal antibody against human MEST (PEG1), suitable for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Useful for studying neurogenesis and adipogenic differentiation.

[Learn More](#)

Parameter	Value
Product Name	MEST Rabbit pAb
Remarks/Alias	PEG1; MEST
Species	Human
GeneID (Human)	4232
GeneID	4232
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 146-335 of human MEST (NP_002393.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5EB52
Protein Weight	39kDa
Shipping	Ice bag