

Anti-Met Polyclonal Antibody For Wb, Ihc-P, Elisa - P08581

Item Number: CM0022844



Introduction

Rabbit polyclonal antibody targeting human MET (Hepatocyte growth factor receptor). Validated for WB, IHC-P, ELISA with cross-reactivity against human, mouse, rat. Suitable for cancer research and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	MET Rabbit pAb
Remarks/Alias	DA11; HGFR; AUTS9; RCCP2; c-Met; DFN897; MET
Species	Human
Gene ID (Human)	4233
Gene ID	4233
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 518-840 of human MET (NP_000236.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08581
Protein Weight	156kDa
Shipping	Ice bag