

Anti-Met Rabbit Polyclonal Antibody For Wb And Elisa - P08581

Item Number: CM0016462



Introduction

CamelBio's Anti-MET Rabbit Polyclonal Antibody (P08581) is validated for Western blotting and ELISA in human and rat samples, targeting the hepatocyte growth factor receptor c-Met for signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MET Rabbit pAb
Remarks/Alias	DA11, HGFR, AUTS9, RCCP2, c-Met, DFNB97, MET
Species	Human
GeneID	4233 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 960-1161 of human MET (NP_000236.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P08581
Protein Weight	156kDa
Shipping	Ice bag