

## Anti-Metap2 Rabbit Polyclonal Antibody For Wb And Elisa - P50579

Item Number: CM0029765



### Introduction

Rabbit polyclonal antibody targeting human METAP2 (methionine aminopeptidase 2), validated for Western blotting and ELISA with cross-reactivity to mouse. Ideal for protein synthesis and translational control studies.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | METAP2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | MAP2; MNPEP; p67eIF2; METAP2                                                                                      |
| Species          | Human                                                                                                             |
| GeneID           | 10988 (Human)                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human METAP2 (NP_006829.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | P50579                                                                                                            |
| Protein Weight   | 53kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |