

Anti-Mettl14 Monoclonal Antibody (Kd Validated) For Wb, Ihc-P, Ip, Elisa - Q9Hce5

Item Number: CM0007364



Introduction

Rabbit monoclonal METTL14 antibody (recombinant immunogen, KD validated) recognizes human/mouse/rat. For WB, IHC-P, IP, ELISA in m6A methylation, circadian clock, and stem cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] METTL14 Rabbit mAb
Remarks/Alias	hMETTL14; [KD Validated] METTL14
Species	Human
GeneID	57721
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human METTL14 (NP_066012.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCE5
Protein Weight	52kDa
Shipping	Ice bag