


Anti-Mettl21C Rabbit Polyclonal Antibody For Wb, Elisa - Q5Vzv1

Item Number: CM0016448



Introduction

High-quality rabbit polyclonal antibody targeting human METTL21C, a lysine methyltransferase. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for methylation and gene regulation studies.

[Learn More](#)

Parameter	Detail
Product Name	METTL21C Rabbit pAb
Remarks/Alias	C13orf39; METTL21C
Species	Human
Gene ID (Human)	196541
Gene ID	196541
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human METTL21C (NP_001010977.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5VZV1
Protein Weight	30kDa
Shipping	Ice bag