

Anti-Mettl2A Polyclonal Antibody For If/Icc, Elisa - Q96Iz6

Item Number: CM0015972

Introduction

Rabbit polyclonal antibody against human METTL2A (tRNA methyltransferase 2A), validated for IF/ICC and ELISA. Cross-reacts with mouse and rat. Suitable for tRNA modification research.

[Learn More](#)

Parameter	Value
Product Name	METTL2A Rabbit pAb
Remarks/Alias	METTL2; METTL2A
Species	Human
GeneID (Human)	339175
GeneID	339175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 139-378 of human METTL2A (NP_859076.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Iz6
Protein Weight	44kDa
Shipping	Ice bag