

Anti-Mettl3 Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ip, If/Icc, Ihc-P, Chip, Elisa - Q86U44

Item Number: CM0009223



Introduction

High-performance recombinant rabbit monoclonal antibody targeting METTL3 (N6-adenosine-methyltransferase catalytic subunit). Validated for WB, IP, IF/ICC, IHC-P, ChIP, and ELISA. Reacts with human, mouse, and rat METTL3. Ideal for epitranscriptomic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] METTL3 Rabbit mAb
Remarks/Alias	M6A; IME4; Spo8; MT-A70; hMETTL3
Species	Human
GeneID (Human)	56339
GeneID	56339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human METTL3 (NP_062826.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86U44
Protein Weight	64 kDa
Shipping	Ice bag