

Anti-Mettl3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q86U44

Item Number: CM0028645



Introduction

Rabbit polyclonal antibody against METTL3 (N6-adenosine-methyltransferase catalytic subunit), validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat METTL3; suitable for RNA methylation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] METTL3 Rabbit pAb
Remarks/Alias	M6A; IME4; Spo8; MT-A70; hMETTL3; METTL3
Species	Human
GeneID (human)	56339
GeneID	56339
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human METTL3 (NP_062826.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86U44
Protein Weight	64kDa
Shipping	Ice bag