

Anti-Mettl4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8N3J2

Item Number: CM0027943

Introduction

METTL4 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets N6-adenine methyltransferase involved in RNA and DNA methylation. Ideal for epigenetic and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	METTL4 Rabbit pAb
Remarks/Alias	HsT661; METTL4
Species	Human
GeneID (Human)	64863
GeneID	64863
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human METTL4 (NP_073751.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N3J2
Protein Weight	54kDa
Shipping	Ice bag