

Anti-Mettl7A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H8H3

Item Number: CM0030669

Introduction

METTL7A rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA applications. Reacts with human, mouse, rat. Targets thiol methyltransferase involved in lncRNA modification and lipid droplet biology. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	METTL7A Rabbit pAb
Remarks/Alias	AAMB; AAM-B; METTL7A
Species	Human
Gene ID (human)	25840
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-244 of human METTL7A (NP_054752.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H8H3
Protein Weight	28kDa
Shipping	Ice bag