

Anti-Mettl7B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Ux53

Item Number: CM0027927

Introduction

Rabbit polyclonal antibody targeting human METTL7B (TMT1B), a thiol S-methyltransferase. Suitable for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for liver metabolism and methyltransferase research.

[Learn More](#)

Parameter	Value
Product Name	METTL7B Rabbit pAb
Remarks/Alias	ALDI; METTL7B
Species	Human
GenelD (Human)	196410
GenelD	196410
Immunogen	Recombinant protein (recombinant fusion protein containing a sequence corresponding to amino acids 165-244 of human METTL7B (NP_689850.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UX53
Protein Weight	28kDa
Shipping	Ice bag