

Anti-Mettl8 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H825

Item Number: CM0022542

Introduction

Rabbit polyclonal antibody specific to human METTL8, validated for WB and ELISA. Reacts with human, mouse, and rat. Suitable for mitochondrial tRNA modification and respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	METTL8 Rabbit pAb
Remarks/Alias	TIP; METTL8
Species	Human
Gene ID (Human)	79828
Gene ID	79828
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human METTL8 (NP_001308089.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H825
Protein Weight	33kDa
Shipping	Ice bag