

Anti-Mettl9 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H1A3

Item Number: CM0020403



Introduction

Rabbit polyclonal antibody against human METTL9, validated for WB and ELISA. Cross-reacts with mouse. Ideal for studying protein histidine methylation and HxH motif targets.

[Learn More](#)

Parameter	Value
Product Name	METTL9 Rabbit pAb
Remarks/Alias	DREV; PAP1; DREV1; CGI-81; hMETTL9; METTL9
Species	Human
GeneID	51108
Immunogen	Recombinant protein of human METTL9
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H1A3
Protein Weight	37kDa
Shipping	Ice bag