

Anti-Mevalonate Kinase (Mvk) Rabbit Monoclonal Antibody For Wb, Elisa - Q03426

Item Number: CM0002180



Introduction

Rabbit monoclonal antibody targeting human mevalonate kinase (MVK, MK, LRBP), validated for WB and ELISA. Suitable for research on cholesterol biosynthesis, sterol metabolism, and related disease pathways.

[Learn More](#)

Parameter	Value
Product Name	MVK Rabbit mAb
Remarks/Alias	MK; LRBP; MVLK; POROK3; MVK
Species	Human
Gene ID (Human)	4598
Gene ID	4598
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 297-396 of human MVK (Q03426)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q03426
Protein Weight	42kDa
Shipping	Ice bag