

Anti-Mex3D Polyclonal Antibody For Wb, Elisa - Q86Xn8

Item Number: CM0021563



Introduction

Rabbit polyclonal antibody targeting human RNA-binding protein MEX3D. Validated for western blot and ELISA, with confirmed cross-reactivity to mouse. Suitable for post-transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	MEX3D Rabbit pAb
Remarks/Alias	MEX3; TINO; RKHD1; MEX-3D; RNF193; OK/SW-cl.4; MEX3D
Species	Human
GeneID (human)	399664
GeneID	399664
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 500-600 of human MEX3D (NP_976049.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q86XN8
Protein Weight	65kDa
Shipping	Ice bag