

Anti-Mfap5 Rabbit Polyclonal Antibody For Wb, Elisa - Q13361

Item Number: CM0022036

Introduction

High-quality rabbit polyclonal antibody targeting human microfibrillar-associated protein 5 (MFAP5), ideal for Western blotting and ELISA applications. Exhibits cross-reactivity with mouse and rat. Supports studies in hematopoiesis, cardiovascular integrity, and elastin-associated microfibrils.

[Learn More](#)

Parameter	Value
Product Name	MFAP5 Rabbit pAb
Remarks/Alias	AAT9; MP25; MAGP2; MAGP-2; MFAP-5; MFAP5
Species	Human
GeneID (Human)	8076
GeneID	8076
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-173 of human MFAP5. (NP_003471.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13361
Protein Weight	20kDa
Shipping	Ice bag