

Anti-Mff Rabbit Polyclonal Antibody (Wb, Ihc-P, Elisa) - Q9Gzy8

Item Number: CM0029138

Introduction

Rabbit polyclonal antibody against human MFF, validated for Western blot, IHC-P, and ELISA. Detects endogenous MFF in human, mouse, and rat samples. Suitable for mitochondrial fission and organelle division studies.

[Learn More](#)

Parameter	Value
Product Name	MFF Rabbit pAb
Remarks/Alias	EMPF2; GL004; C2orf33; MFF
Species	Human
GeneID (Human)	56947
GeneID	56947
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MFF (NP_064579.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZY8
Protein Weight	38kDa
Shipping	Ice bag