

Anti-Mfg-E8 Monoclonal Antibody For Wb, If/Icc, Elisa - Q08431

Item Number: CM0010273

Introduction

MFG-E8 Rabbit mAb validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets lactadherin (Q08431), a milk fat globule-EGF factor 8 involved in apoptotic cell clearance, angiogenesis, and mucosal healing. Ideal for cancer and fertilization research.

[Learn More](#)

Parameter	Value
Product Name	MFG-E8 Rabbit mAb
Remarks/Alias	BA46; HMFG; MFGM; SED1; hP47; EDIL1; MFG-E8; SPAG10; OAcGD3S; HsT19888
Species	human
GeneID (human)	4240
GeneID	4240
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 231-387 of human MFG-E8 (NP_005919.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08431
Protein Weight	43kDa
Shipping	Ice bag