

## Anti-Mfge8 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q08431

Item Number: CM0017758

### Introduction

MFGE8 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IHC-P, and ELISA, cross-reacts with human, mouse, and rat. For apoptosis, angiogenesis, mucosal healing, and gamete interaction research. Immunogen: aa 1-200; ~43 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | MFGE8 Rabbit pAb                                                                                                 |
| Remarks/Alias    | BA46; HMFG; MFGM; SED1; hP47; EDIL1; MFG-E8; SPAG10; OAcGD3S; HsT19888; MFGE8                                    |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 4240                                                                                                             |
| Gene ID          | 4240                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human MFGE8 (NP_005919.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q08431                                                                                                           |
| Protein Weight   | 43 kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                          |