

Anti-Mfi2 (Melanotransferrin) Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P08582

Item Number: CM0004549



Introduction

Rabbit monoclonal antibody targeting human melanotransferrin (MFI2/CD228), suitable for WB, FC, and ELISA. Cross-reacts with rat. Ideal for iron metabolism and melanoma research.

[Learn More](#)

Parameter	Value
Product Name	MFI2 Rabbit mAb
Remarks/Alias	MTf; MFI2; MTF1; CD228; MAP97
Species	Human
GeneID (Human)	4241
GeneID	4241
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-708 of human MFI2 (NP_005920.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human, Rat
SWISS	P08582
Protein Weight	32kDa/80kDa
Shipping	Ice bag