

Anti-Mfi2 (Melanotransferrin) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08582

Item Number: CM0027107



Introduction

Rabbit polyclonal antibody against human MFI2 (melanotransferrin, CD228). Validated for Western blot, immunofluorescence, and ELISA with cross-reactivity to mouse and rat. Suitable for iron transport and melanoma research.

[Learn More](#)

Parameter	Value
Product Name	MFI2 Rabbit pAb
Remarks/Alias	MTf; MFI2; MTF1; CD228; MAP97
Species	Human
GeneID (Human)	4241
GeneID	4241
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-200 of human MFI2 (NP_201573.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08582
Protein Weight	80kDa
Shipping	Ice bag