


Anti-Mfsd2A Polyclonal Antibody For Wb, Elisa - Q8Na29

Item Number: CM0020407



Introduction

Rabbit polyclonal antibody against human MFSD2A, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for BBB and lipid transport research.

[Learn More](#)

Parameter	Value
Product Name	MFSD2A Rabbit pAb
Remarks/Alias	NLS1; MFSD2; MCPH15; SLC59A1; HsMFSD2A; NEDMISBA; MFSD2A
Species	Human
GeneID (human)	84879
GeneID	84879
Immunogen	Recombinant fusion protein containing a sequence within amino acids 464-543 of human MFSD2A (NP_001129965.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	Q8NA29
Protein Weight	60 kDa
Shipping	Ice bag