


Anti-Mfsd2A/Nls1 Rabbit Monoclonal Antibody For Wb, Elisa - Q8Na29

Item Number: CM0006182



Introduction

Rabbit monoclonal antibody targeting human MFSD2A/NLS1 (Sodium-dependent lysophosphatidylcholine symporter 1). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for blood-brain barrier and lipid transport studies.

[Learn More](#)

Parameter	Value
Product Name	MFSD2A/NLS1 Rabbit mAb
Remarks/Alias	NLS1; MFSD2; MCPH15; SLC59A1; HsMFSD2A; NEDMISBA; MFSD2A/NLS1
Species	Human
GeneID	84879
Immunogen	Recombinant protein of human MFSD2A/NLS1
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NA29
Protein Weight	60kDa
Shipping	Ice bag