

Anti-Mgat1 Polyclonal Antibody For Wb, Elisa - P26572

Item Number: CM0014554

Introduction

High-quality rabbit polyclonal antibody targeting MGAT1 (GNT-I), ideal for WB and ELISA applications. Shows cross-reactivity with human and rat samples. Essential tool for studying N-glycan biosynthesis and Golgi function.

[Learn More](#)

Parameter	Value
Product Name	MGAT1 Rabbit pAb
Remarks/Alias	GnTI; MGAT; GLCT1; GLYT1; GNT-1; GNT-I; GLCNAC-TI; MGAT1
Species	Human
GeneID (Human)	4245
GeneID	4245
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MGAT1 (NP_002397.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P26572
Protein Weight	51kDa
Shipping	Ice bag