

Anti-Mgat2 Rabbit Polyclonal Antibody For Wb And Elisa - Q10469

Item Number: CM0026843



Introduction

Rabbit polyclonal antibody targeting human MGAT2 (GlcNAc-T II), essential for protein N-glycosylation, validated in WB and ELISA across human, mouse, and rat samples. Ideal for glycobiology research.

[Learn More](#)

Parameter	Value
Product Name	MGAT2 Rabbit pAb
Remarks/Alias	GNT2; CDG2A; CDGS2; GNT-II; GLCNACTII; MGAT2
Species	Human
Gene ID (Human)	4247
Gene ID	4247
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 178-447 of human MGAT2 (NP_002399.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q10469
Protein Weight	52kDa
Shipping	Ice bag