

Anti-Mgat3 Monoclonal Antibody For Wb, Elisa - Q09327

Item Number: CM0006714

Introduction

Rabbit monoclonal antibody against human MGAT3 (GNT-III), validated for Western Blot and ELISA. Cross-reacts with mouse. Suitable for glycoprotein biosynthesis and N-glycan branching studies.

[Learn More](#)

Parameter	Value
Product Name	MGAT3 Rabbit mAb
Remarks/Alias	GNT3; GNT-III; MGAT3
Species	Human
Gene ID	4248
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 381-480 of human MGAT3 (NP_002400.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q09327
Protein Weight	61kDa
Shipping	Ice bag