

Anti-Mgat3 Polyclonal Antibody For Wb And Elisa - Q09327

Item Number: CM0027264

Introduction

Rabbit polyclonal antibody against human MGAT3 (GNT-III) validated for Western blot and ELISA. Detects MGAT3 in human, mouse, and rat samples; involved in N-glycan biosynthesis and cell adhesion regulation.

[Learn More](#)

Parameter	Value
Product Name	MGAT3 Rabbit pAb
Remarks/Alias	GNT3; GNT-III; MGAT3
Species	Human
Gene ID (Human)	4248
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 364-533 of human MGAT3 (NP_002400.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q09327
Protein Weight	61kDa
Shipping	Ice bag