


Anti-Mgat4A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Um21

Item Number: CM0017852



Introduction

Rabbit polyclonal antibody targeting human MGAT4A (GNT-IVA), validated for Western blot and ELISA. Ideal for studying glycosylation, glucose transport, and pancreatic beta cell function. Cross-reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	MGAT4A Rabbit pAb
Remarks/Alias	GNT-IV; GnT-4a; GNT-IVA; MGAT4A
Species	Human
GeneID (Human)	11320
GeneID	11320
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-100 of human MGAT4A (NP_036346.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UM21
Protein Weight	62kDa
Shipping	Ice bag