

Anti-Mgea5 Polyclonal Antibody For Wb, Ihc-P, Elisa - O60502

Item Number: CM0029644



Introduction

Rabbit polyclonal antibody targeting human MGEA5/OGA. Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Suitable for studying O-GlcNAcylation and related diseases. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	MGEA5 Rabbit pAb
Remarks/Alias	MEA5; MGEA5; NCOAT
Species	Human
Gene ID (Human)	10724
Gene ID	10724
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-580 of human MGEA5 (NP_001135906.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60502
Protein Weight	103kDa
Shipping	Ice bag