

Anti-Mgea5/Oga Monoclonal Antibody For Wb, Ihc-P, Elisa - O60502

Item Number: CM0005046



Introduction

Anti-MGEA5/OGA monoclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for O-GlcNAcase research and IVD assay development. Provided by CamelBio's comprehensive IVD raw material supply platform.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MGEA5/OGA Rabbit mAb
Remarks/Alias	MEA5; MGEA5; NCOAT
Species	Human
Gene ID (Human)	10724
Gene ID	10724
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 433-633 of human MGEA5/OGA NP_036347.1
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60502
Protein Weight	103 kDa
Shipping	Ice bag